


# PARI PEP®-Systems

## Time-saving combination of physiotherapy and inhalation treatment

### The key features at a glance

- Mobilises secretions and thus facilitates clearing of the lungs
- Stabilises the respiratory tract by adjustable expiration resistance: counteracts bronchial collapse
- Pressure display for setting expiration resistance and for monitoring during PEP treatment (PARI PEP® System II)
- Can also be used without the nebuliser in physiotherapy



① **PARI PEP® S System** for **PARI** Nebuliser (except PARI LC SPRINT BABY Nebuliser, PARI BABY Nebuliser and PARI TIA Nebuliser S)  
Item No. 018G4000

② **PARI PEP® System I** for **PARI LC PLUS®**, **PARI LC STAR®** and **PARI LL Nebuliser**.  
Without manometer  
Item No. 018G6101

② **PARI PEP® System II** for **PARI LC PLUS®**, **PARI LC STAR®** and **PARI LL Nebuliser**.  
With manometer  
Item No. 018G6200

### Great hygiene benefit - the PARI PEP® Systems can be:

- washed in a dishwasher
- disinfected by boiling
- chemically disinfected
- autoclaved

### Accessories for the PARI PEP® Systems



**Manometer**  
0-100 mbar, including Pressure Hose  
Item No. 018B1600

**PARI PEP® Tubing adapter**  
Pack of 3, made of silicon, with one adapter for manometer  
Item No. 041G4580

